

Nailing System



MEDI ORTHO LINK

Bringing Life Back To Your Bones



Interlocking Nail For Tibia



Tibia Nail (Gundrill)

Cat No. 155.13 / T155.13

Stainless Steel

Cat No.	Dia.	Length	Cat No.	Dia.	Length
155.13.0828	08 MM	28 CM	155.13.1028	10 MM	28 CM
155.13.0829	08 MM	29 CM	155.13.1029	10 MM	29 CM
155.13.0830	08 MM	30 CM	155.13.1030	10 MM	30 CM
155.13.0831	08 MM	31 CM	155.13.1031	10 MM	31 CM
155.13.0832	08 MM	32 CM	155.13.1032	10 MM	32 CM
155.13.0833	08 MM	33 CM	155.13.1033	10 MM	33 CM
155.13.0834	08 MM	34 CM	155.13.1034	10 MM	34 CM
155.13.0835	08 MM	35 CM	155.13.1035	10 MM	35 CM
155.13.0836	08 MM	36 CM	155.13.1036	10 MM	36 CM
155.13.0837	08 MM	37 CM	155.13.1037	10 MM	37 CM
155.13.0838	08 MM	38 CM	155.13.1038	10 MM	38 CM
155.13.0928	09 MM	28 CM	155.13.1128	11 MM	28 CM
155.13.0929	09 MM	29 CM	155.13.1129	11 MM	29 CM
155.13.0930	09 MM	30 CM	155.13.1130	11 MM	30 CM
155.13.0931	09 MM	31 CM	155.13.1131	11 MM	31 CM
155.13.0932	09 MM	32 CM	155.13.1132	11 MM	32 CM
155.13.0933	09 MM	33 CM	155.13.1133	11 MM	33 CM
155.13.0934	09 MM	34 CM	155.13.1134	11 MM	34 CM
155.13.0935	09 MM	35 CM	155.13.1135	11 MM	35 CM
155.13.0936	09 MM	36 CM	155.13.1136	11 MM	36 CM
155.13.0937	09 MM	37 CM	155.13.1137	11 MM	37 CM
155.13.0938	09 MM	38 CM	155.13.1138	11 MM	38 CM

Titanium

Cat No.	Dia.	Length
T155.0828	08 MM	28 CM
T155.0830	08 MM	30 CM
T155.0832	08 MM	32 CM
T155.0834	08 MM	34 CM
T155.0836	08 MM	36 CM
T155.0838	08 MM	38 CM
T155.0928	09 MM	28 CM
T155.0930	09 MM	30 CM
T155.0932	09 MM	32 CM
T155.0934	09 MM	34 CM
T155.0936	09 MM	36 CM
T155.0938	09 MM	38 CM
T155.1028	10 MM	28 CM
T155.1030	10 MM	30 CM
T155.1032	10 MM	32 CM
T155.1034	10 MM	34 CM
T155.1036	10 MM	36 CM
T155.1038	10 MM	38 CM



Interlocking Nail For Femur



Femur Nail (Gundrill)

Cat No: 156.13 / T156.13

Stainless Steel

Cat No.	Dia.	Length	Cat No.	Dia.	Length
156.13.0934	09 MM	34 CM	156.13.1134	11 MM	34 CM
156.13.0936	09 MM	36 CM	156.13.1136	11 MM	36 CM
156.13.0938	09 MM	38 CM	156.13.1138	11 MM	38 CM
156.13.0940	09 MM	40 CM	156.13.1140	11 MM	40 CM
156.13.0942	09 MM	42 CM	156.13.1142	11 MM	42 CM
156.13.0944	09 MM	44 CM	156.13.1144	11 MM	44 CM
156.13.1034	10 MM	34 CM	156.13.1234	12 MM	34 CM
156.13.1036	10 MM	36 CM	156.13.1236	12 MM	36 CM
156.13.1038	10 MM	38 CM	156.13.1238	12 MM	38 CM
156.13.1040	10 MM	40 CM	156.13.1240	12 MM	40 CM
156.13.1042	10 MM	42 CM	156.13.1242	12 MM	42 CM
156.13.1044	10 MM	44 CM	156.13.1244	12 MM	44 CM

Titanium

Cat No.	Dia.	Length
T156.0934	09 MM	34 CM
T156.0936	09 MM	36 CM
T156.0938	09 MM	38 CM
T156.0940	09 MM	40 CM
T156.0942	09 MM	42 CM
T156.0944	09 MM	44 CM
T156.1034	10 MM	34 CM
T156.1036	10 MM	36 CM
T156.1038	10 MM	38 CM
T156.1040	10 MM	40 CM
T156.1042	10 MM	42 CM
T156.1044	10 MM	44 CM
T156.1134	11 MM	34 CM
T156.1136	11 MM	36 CM
T156.1138	11 MM	38 CM
T156.1140	11 MM	40 CM
T156.1142	11 MM	42 CM
T156.1144	11 MM	44 CM



Cat. No. 150, T150

I.L. Bolt 3.9 mm*



Thread Dia. : 3.9 MM
 Core Dia. : 3.2 MM
 Head Dia. : 8.0 MM
 Drill Bit Use : 3.2 MM
 Screw Driver Use : 4.5 MM

Cat. No. 151, T151

I.L. Bolt 4.9 mm



Thread Dia. : 4.9 MM
 Core Dia. : 4.3 MM
 Head Dia. : 8.0 MM
 Drill Bit Use : 4.0 MM
 Screw Driver Use : 4.5 MM

Cat. No. 183 / T183 Cap For Tibia & Femur Nail

Size : 00 MM, 05 MM



* Use Only 08 mm Distal Hole.

Instrument Set For Tibia & Femur Interlocking Nail



Aiming Arm For Tibia & Femur Nail



Jig For Tibia & Femur Nail



Conical Bolt



Protection Sleeve



Pointer



Drill Sleeve for 4.0 MM Drill Bit



Drill Sleeve for 3.2 MM Drill Bit



T - Handle (Wrench)



Spanner (Ordinary)
12 X 13, 18 X 13



Extractor Rod



Round Hammer



Impactor / Extractor Handle



Depth Gauge Long



Bone Awl



Tissue Protector



Proximal Reamer 13 MM



Screw Driver 4.5 X 200 MM



Drill Bits 4.0 MM X 10°



Drill Bits 3.2 MM X 08°

Interlocking Nail For Tibia With Distal Jig



Tibia Nail With Distal Jig

Cat No.:- 155.D / T155.D

Stainless Steel

Cat No.	Dia.	Length	Cat No.	Dia.	Length
155.D.08285	08 MM	285 MM	155.D.10285	10 MM	285 MM
155.D.08300	08 MM	300 MM	155.D.10300	10 MM	300 MM
155.D.08315	08 MM	315 MM	155.D.10315	10 MM	315 MM
155.D.08330	08 MM	330 MM	155.D.10330	10 MM	330 MM
155.D.08345	08 MM	345 MM	155.D.10345	10 MM	345 MM
155.D.08360	08 MM	360 MM	155.D.10360	10 MM	360 MM
155.D.08375	08 MM	375 MM	155.D.10375	10 MM	375 MM
155.D.08390	08 MM	390 MM	155.D.10390	10 MM	390 MM
155.D.08405	08 MM	405 MM	155.D.10405	10 MM	405 MM
155.D.09285	09 MM	285 MM	155.D.11285	11 MM	285 MM
155.D.09300	09 MM	300 MM	155.D.11300	11 MM	300 MM
155.D.09315	09 MM	315 MM	155.D.11315	11 MM	315 MM
155.D.09330	09 MM	330 MM	155.D.11330	11 MM	330 MM
155.D.09345	09 MM	345 MM	155.D.11345	11 MM	345 MM
155.D.09360	09 MM	360 MM	155.D.11360	11 MM	360 MM
155.D.09375	09 MM	375 MM	155.D.11375	11 MM	375 MM
155.D.09390	09 MM	390 MM	155.D.11390	11 MM	390 MM
155.D.09405	09 MM	405 MM	155.D.11405	11 MM	405 MM

Titanium

Cat No.	Dia.	Length
T155.D.08285	08 MM	285 MM
T155.D.08300	08 MM	300 MM
T155.D.08315	08 MM	315 MM
T155.D.08330	08 MM	330 MM
T155.D.08345	08 MM	345 MM
T155.D.08360	08 MM	360 MM
T155.D.08375	08 MM	375 MM
T155.D.08390	08 MM	390 MM
T155.D.08405	08 MM	405 MM
T155.D.09285	09 MM	285 MM
T155.D.09300	09 MM	300 MM
T155.D.09315	09 MM	315 MM
T155.D.09330	09 MM	330 MM
T155.D.09345	09 MM	345 MM
T155.D.09360	09 MM	360 MM
T155.D.09375	09 MM	375 MM
T155.D.09390	09 MM	390 MM
T155.D.09405	09 MM	405 MM
T155.D.10285	10 MM	285 MM
T155.D.10300	10 MM	300 MM
T155.D.10315	10 MM	315 MM
T155.D.10330	10 MM	330 MM
T155.D.10345	10 MM	345 MM
T155.D.10360	10 MM	360 MM
T155.D.10375	10 MM	375 MM
T155.D.10390	10 MM	390 MM
T155.D.10405	10 MM	405 MM

Cat. No. 150 / T150

I.L. Bolt 3.9 MM*



Product Profile

Thread Dia.	: 3.9 MM
Core Dia.	: 3.2 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 3.2 MM
Screw Driver Use	: 4.5 MM

Cat. No. 151 / T151

I.L. Bolt 4.9 MM



Product Profile

Thread Dia.	: 4.9 MM
Core Dia.	: 4.3 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 4.0 MM
Screw Driver Use	: 4.5 MM



Interlocking Nail For Femur With Distal Jig



Femur Nail With Distal Jig

Cat No. :- 156.D / T156.D

Stainless Steel

Cat No.	Dia.	Length	Cat No.	Dia.	Length
156.D.09330	09 MM	330 MM	156.D.11330	11 MM	330 MM
156.D.09345	09 MM	345 MM	156.D.11345	11 MM	345 MM
156.D.09360	09 MM	360 MM	156.D.11360	11 MM	360 MM
156.D.09375	09 MM	375 MM	156.D.11375	11 MM	375 MM
156.D.09390	09 MM	390 MM	156.D.11390	11 MM	390 MM
156.D.09405	09 MM	405 MM	156.D.11405	11 MM	405 MM
156.D.09420	09 MM	420 MM	156.D.11420	11 MM	420 MM
156.D.09435	09 MM	435 MM	156.D.11435	11 MM	435 MM
156.D.09450	09 MM	450 MM	156.D.11450	11 MM	450 MM
156.D.10330	10 MM	330 MM	156.D.12330	12 MM	330 MM
156.D.10345	10 MM	345 MM	156.D.12345	12 MM	345 MM
156.D.10360	10 MM	360 MM	156.D.12360	12 MM	360 MM
156.D.10375	10 MM	375 MM	156.D.12375	12 MM	375 MM
156.D.10390	10 MM	390 MM	156.D.12390	12 MM	390 MM
156.D.10405	10 MM	405 MM	156.D.12405	12 MM	405 MM
156.D.10420	10 MM	420 MM	156.D.12420	12 MM	420 MM
156.D.10435	10 MM	435 MM	156.D.12435	12 MM	435 MM
156.D.10450	10 MM	450 MM	156.D.12450	12 MM	450 MM

Titanium

Cat No.	Dia.	Length
T156.D.09330	09 MM	330 MM
T156.D.09345	09 MM	345 MM
T156.D.09360	09 MM	360 MM
T156.D.09375	09 MM	375 MM
T156.D.09390	09 MM	390 MM
T156.D.09405	09 MM	405 MM
T156.D.09420	09 MM	420 MM
T156.D.09435	09 MM	435 MM
T156.D.09450	09 MM	450 MM
T156.D.10330	10 MM	330 MM
T156.D.10345	10 MM	345 MM
T156.D.10360	10 MM	360 MM
T156.D.10375	10 MM	375 MM
T156.D.10390	10 MM	390 MM
T156.D.10405	10 MM	405 MM
T156.D.10420	10 MM	420 MM
T156.D.10435	10 MM	435 MM
T156.D.10450	10 MM	450 MM
T156.D.11330	11 MM	330 MM
T156.D.11345	11 MM	345 MM
T156.D.11360	11 MM	360 MM
T156.D.11375	11 MM	375 MM
T156.D.11390	11 MM	390 MM
T156.D.11405	11 MM	405 MM
T156.D.11420	11 MM	420 MM
T156.D.11435	11 MM	435 MM
T156.D.11450	11 MM	450 MM



Cat. No. 183 / T183
Cap For Tibia & Femur Nail With Distal Jig
Size : 00 MM, 05 MM



Instrument Set For Tibia & Femur Interlocking Nail With Distal Jig

Aiming Arm For Tibia & Femur Nail



Distal Jig For Tibia & Femur Nail



Conical Bolt



Protection Sleeve For Tibia Jig



Protection Sleeve For Femur Jig



Pointer



Drill Sleeve for 4.0 MM Drill Bit



Drill Sleeve for 3.2 MM Drill Bit



T - Handle (Wrench)



Spanner (Ordinary)
12 X 13, 18 X 13



Extractor Rod



Round Hammer



Impactor / Extractor Handle



Depth Gauge Long



Bone Awl



Tissue Protector



Proximal Reamer 13 MM



Screw Driver 4.5 X 200 MM



Drill Bits 4.0 MM X 10"



Drill Bits 3.2 MM X 10"



Interlocking Nail For Humerus Solid (With Compression System)



Humerus Nail Solid
Cat. No. 154 / T154

Stainless Steel

Cat No.	Dia.	Length
154.620	06 MM	20 CM
154.621	06 MM	21 CM
154.622	06 MM	22 CM
154.623	06 MM	23 CM
154.624	06 MM	24 CM
154.625	06 MM	25 CM
154.626	06 MM	26 CM
154.627	06 MM	27 CM
154.628	06 MM	28 CM
154.629	06 MM	29 CM
154.630	06 MM	30 CM

Stainless Steel

Cat No.	Dia.	Length
154.720	07 MM	20 CM
154.721	07 MM	21 CM
154.722	07 MM	22 CM
154.723	07 MM	23 CM
154.724	07 MM	24 CM
154.725	07 MM	25 CM
154.726	07 MM	26 CM
154.727	07 MM	27 CM
154.728	07 MM	28 CM
154.729	07 MM	29 CM
154.730	07 MM	30 CM
154.820	08 MM	20 CM
154.821	08 MM	21 CM
154.822	08 MM	22 CM
154.823	08 MM	23 CM
154.824	08 MM	24 CM
154.825	08 MM	25 CM
154.826	08 MM	26 CM
154.827	08 MM	27 CM
154.828	08 MM	28 CM
154.829	08 MM	29 CM
154.830	08 MM	30 CM

Titanium

Cat No.	Dia.	Length
T154.620	06 MM	20 CM
T154.621	06 MM	21 CM
T154.622	06 MM	22 CM
T154.623	06 MM	23 CM
T154.624	06 MM	24 CM
T154.625	06 MM	25 CM
T154.626	06 MM	26 CM
T154.627	06 MM	27 CM
T154.628	06 MM	28 CM
T154.629	06 MM	29 CM
T154.630	06 MM	30 CM



Cat. No. 153 / T153
Cap For Humerus Nail

Size : 00 MM, 05 MM, 10 MM, 15 MM



Interlocking Nail For Humerus Cannulated (With Compression System)



Humerus Nail Cannulated Cat. No. 154.C / T154.C

Stainless Steel

Cat No.	Dia.	Length
154C.720	07 MM	20 CM
154C.721	07 MM	21 CM
154C.722	07 MM	22 CM
154C.723	07 MM	23 CM
154C.724	07 MM	24 CM
154C.725	07 MM	25 CM
154C.726	07 MM	26 CM
154C.727	07 MM	27 CM
154.C728	07 MM	28 CM
154C.729	07 MM	29 CM
154C.730	07 MM	30 CM
154C.820	08 MM	20 CM
154C.821	08 MM	21 CM
154C.822	08 MM	22 CM
154C.823	08 MM	23 CM
154C.824	08 MM	24 CM
154C.825	08 MM	25 CM
154C.826	08 MM	26 CM
154C.827	08 MM	27 CM
154C.828	08 MM	28 CM
154C.829	08 MM	29 CM
154C.830	08 MM	30 CM

Titanium

Cat No.	Dia.	Length
T154C.720	07 MM	20 CM
T154C.721	07 MM	21 CM
T154C.722	07 MM	22 CM
T154C.723	07 MM	23 CM
T154C.724	07 MM	24 CM
T154C.725	07 MM	25 CM
T154C.726	07 MM	26 CM
T154C.727	07 MM	27 CM
T154.C728	07 MM	28 CM
T154C.729	07 MM	29 CM
T154C.730	07 MM	30 CM
T154C.820	08 MM	20 CM
T154C.821	08 MM	21 CM
T154C.822	08 MM	22 CM
T154C.823	08 MM	23 CM
T154C.824	08 MM	24 CM
T154C.825	08 MM	25 CM
T154C.826	08 MM	26 CM
T154C.827	08 MM	27 CM
T154C.828	08 MM	28 CM
T154C.829	08 MM	29 CM
T154C.830	08 MM	30 CM



Cat. No. 148 / T148

TIT I.L. Bolt 2.9 MM*



Product Profile

Thread Dia.	: 2.9 MM
Core Dia.	: 2.4 MM
Head Dia.	: 6.0 MM
Drill Bit Use	: 2.5 MM
Screw Driver	
Use	: 3.5 MM

Cat. No. 149 / T149

TIT I.L. Bolt 3.4 MM



Product Profile

Thread Dia.	: 3.4 MM
Core Dia.	: 3.2 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 2.7 MM
Screw Driver	
Use	: 3.5 MM

Cat. No. 150 / T150

TIT I.L. Bolt 3.9 MM*



Product Profile

Thread Dia.	: 3.9 MM
Core Dia.	: 3.2 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 3.2 MM
Screw Driver	
Use	: 4.5 MM

* Use Only 06 mm Distal Hole.

* Use Only Dynamic Hole.

Instrument Set For Humerus Interlocking Nail



Aiming Arm For Humerus Nail



Jig For Humerus Nail



Conical Bolt



Sleeve for Extraction



Compression Device



Protection Sleeve



Pointer



Drill Sleeve for 3.2 MM Drill Bit



Drill Sleeve for 2.7 MM Drill Bit



Spanner (Ordinary)
12 X 13

Extractor Rod



Round Hammer



Small Hammer Unit



Depth Gauge Long



Cannulated Bone Awl



Tissue Protector



Proximal Reamer For Humerus



Screw Driver 3.5 X 200 MM



Screw Driver 4.5 X 200 MM



Drill Bits 3.2 MM X 08"



Drill Bits 2.7 MM X 08"



Interlocking Nail For Proximal Femoral 130° & 135° (P.F.N.)



Proximal Femoral Nail (P.F.N.)
(L & R) (130° & 135°)
Cat No.:- 180.L/R / T180.L/R

Proximal Femoral Nail (P.F.N.)
Short (130° & 135°)
Cat No.:- 180 / T180

Stainless Steel

Cat No.	Dia	Length	Cat No.	Dia	Length
180.L/R.0934	09 MM	34 CM	180.L/R.1134	11 MM	34 CM
180.L/R.0936	09 MM	36 CM	180.L/R.1136	11 MM	36 CM
180.L/R.0938	09 MM	38 CM	180.L/R.1138	11 MM	38 CM
180.L/R.0940	09 MM	40 CM	180.L/R.1140	11 MM	40 CM
180.L/R.0942	09 MM	42 CM	180.L/R.1142	11 MM	42 CM
180.L/R.0944	09 MM	44 CM	180.L/R.1144	11 MM	44 CM
180.L/R.1034	10 MM	34 CM	180.L/R.1234	12 MM	34 CM
180.L/R.1036	10 MM	36 CM	180.L/R.1236	12 MM	36 CM
180.L/R.1038	10 MM	38 CM	180.L/R.1238	12 MM	38 CM
180.L/R.1040	10 MM	40 CM	180.L/R.1240	12 MM	40 CM
180.L/R.1042	10 MM	42 CM	180.L/R.1242	12 MM	42 CM
180.L/R.1044	10 MM	44 CM	180.L/R.1244	12 MM	44 CM

Titanium

Cat No.	Dia	Length	Cat No.	Dia	Length
T180.L/R.0934	09 MM	34 CM	T180.L/R.1134	11 MM	34 CM
T180.L/R.0936	09 MM	36 CM	T180.L/R.1136	11 MM	36 CM
T180.L/R.0938	09 MM	38 CM	T180.L/R.1138	11 MM	38 CM
T180.L/R.0940	09 MM	40 CM	T180.L/R.1140	11 MM	40 CM
T180.L/R.0942	09 MM	42 CM	T180.L/R.1142	11 MM	42 CM
T180.L/R.0944	09 MM	44 CM	T180.L/R.1144	11 MM	44 CM
T180.L/R.1034	10 MM	34 CM	T180.L/R.1234	12 MM	34 CM
T180.L/R.1036	10 MM	36 CM	T180.L/R.1236	12 MM	36 CM
T180.L/R.1038	10 MM	38 CM	T180.L/R.1238	12 MM	38 CM
T180.L/R.1040	10 MM	40 CM	T180.L/R.1240	12 MM	40 CM
T180.L/R.1042	10 MM	42 CM	T180.L/R.1242	12 MM	42 CM
T180.L/R.1044	10 MM	44 CM	T180.L/R.1244	12 MM	44 CM

Stainless Steel

Cat No.	Dia.	Length
180.0925	09 MM	25 CM
180.1025	10 MM	25 CM
180.1125	11 MM	25 CM
180.1225	12 MM	25 CM

Titanium

Cat No.	Dia.	Length.
T180.0925	09 MM	25 CM
T180.1025	10 MM	25 CM
T180.1125	11 MM	25 CM
T180.1225	12 MM	25 CM



Cat. No. 182.P / T182.P Cap For P.F.N. & T.F.N. Nail
Size : 00 MM, 05 MM, 10 MM



Interlocking Nail For Trochanteric Femoral Nail 130° & 135° (T.F.N.)



Trochanteric Femoral Nail
130° & 135° (T.F.N.)
Cat No. :- 189 / T189

Stainless Steel

Cat No.	Dia.	Length
189.0918	09 MM	18 CM
189.1018	10 MM	18 CM
189.1118	11 MM	18 CM
189.1218	12 MM </td <td>18 CM</td>	18 CM

Titanium

Cat No.	Dia.	Length
T 189.0918	09 MM	18 CM
T 189.1018	10 MM	18 CM
T 189.1118	11 MM	18 CM
T 189.1218	12 MM	18 CM



Cat. No. 151 / T151
TIT I.L. Bolt 4.9 MM



Thread Dia. : 4.9 MM
Core Dia. : 4.3 MM
Head Dia. : 8.0 MM
Drill Bit Use : 4.0 MM
Screw Driver Use : 4.5 MM

Cat. No. 178 / T178
Proximal Bolt 8.0 MM



Product Profile

Thread Dia. : 8.0 MM
Core Dia. : 7.9 MM
Head Dia. : 10.0 MM
Screw Driver Use : 5.0 MM

Cat. No. 179 / T179
Proximal Bolt 6.4 MM



Product Profile

Thread Dia. : 6.4 MM
Core Dia. : 4.8 MM
Head Dia. : 8.0 MM
Screw Driver Use : 5.0 MM

Instrument Set For P.F.N. / T.F.N. Interlocking Nail



Aiming Arm For P.F.N / T.F.N Nail 130°
Aiming Arm For P.F.N / T.F.N Nail 135°



Combine Jig For P.F.N / T.F.N Nail 130° & 135°

Conical Bolt



Protection Sleeve For 8.0 MM Screw



Protection Sleeve (2 Nos.) For 6.4 MM Screw



Drill Sleeve for 4.0 MM Drill Bit



Guide Wire Sleeve for 2.0 MM
Use For 6.4 MM Bolt



Guide Wire Sleeve for 2.0 MM
Use For 8.0 MM Bolt



Pointer Long



T - Handle (Wrench)



Spanner (Ordinary)
12 X 13 , 18 X 13



Extractor Rod



Round Hammer



Impactor / Extractor Handle



Introducer Hammer Unit



Depth Gauge Extra Long



Bone Awl



Tissue Protector



Proximal Reamer 15 MM



Screw Driver 4.5 X 250 MM



C. C. Screw Driver 5.0 X 250 MM



C.C.Drill Bit Q/C 8.0 MM



C.C.Drill Bit Q/C 6.4 MM




Drill Bits 4.0 MM X 10"



T Handle For Drill Bits





Interlocking
Multifunction Tibia Nail
(MTN)

Interlocking Nail For Multifunction Tibia Nail (MTN)



Multifunction Tibia Nail (MTN)

Cat No.- 186 / T186

Stainless Steel

Cat No.	Dia.	Length
186.0828	08 MM	28 CM
186.0830	08 MM	30 CM
186.0832	08 MM	32 CM
186.0834	08 MM	34 CM
186.0836	08 MM	36 CM
186.0838	08 MM	38 CM
186.0928	09 MM	28 CM
186.0930	09 MM	30 CM
186.0932	09 MM	32 CM
186.0934	09 MM	34 CM
186.0936	09 MM	36 CM
186.0938	09 MM	38 CM
186.1028	10 MM	28 CM
186.1030	10 MM	30 CM
186.1032	10 MM	32 CM
186.1034	10 MM	34 CM
186.1036	10 MM	36 CM
186.1038	10 MM	38 CM

Titanium

Cat No.	Dia.	Length
T186.0828	08 MM	28 CM
T186.0830	08 MM	30 CM
T186.0832	08 MM	32 CM
T186.0834	08 MM	34 CM
T186.0836	08 MM	36 CM
T186.0838	08 MM	38 CM
T186.0928	09 MM	28 CM
T186.0930	09 MM	30 CM
T186.0932	09 MM	32 CM
T186.0934	09 MM	34 CM
T186.0936	09 MM	36 CM
T186.0938	09 MM	38 CM
T186.1028	10 MM	28 CM
T186.1030	10 MM	30 CM
T186.1032	10 MM	32 CM
T186.1034	10 MM	34 CM
T186.1036	10 MM	36 CM
T186.1038	10 MM	38 CM



Cat. No. 150 / T150

I.L. Bolt 3.9 MM *



Product Profile

Thread Dia.	: 3.9 MM
Core Dia.	: 3.2 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 3.2 MM
Screw Driver Use	: 4.5 MM

Cat. No. 151 / T151

I.L. Bolt 4.9 MM



Product Profile

Thread Dia.	: 4.9 MM
Core Dia.	: 4.3 MM
Head Dia.	: 8.0 MM
Drill Bit Use	: 4.0 MM
Screw Driver Use	: 4.5 MM

Cat. No. 182.M / T182.M Cap For MTN

Size : 00 MM, 05 MM



* Use Only 08 & 09 mm Distal Hole.

Instrument Set For Multifunction Tibia Nail (MTN)



Aiming Arm For Multifunction Tibia Nail



Jig For Multifunction Tibia Nail



Conical Bolt



Protection Sleeve



Pointer



Drill Sleeve for 4.0 MM Drill Bit



T - Handle (Wrench)



Spanner (Ordinary)
12 X 13 , 18 X 13



Extractor Rod

Round Hammer



Impactor / Extractor Handle



Depth Gauge Long



Bone Awl



Tissue Protector



Proximal Reamer 11 MM



Screw Driver 4.5 X 200 MM



Drill Bits 4.0 MM X 10"



Drill Bits 3.2 MM X 08"



A photograph of a man and a woman walking barefoot on a sandy beach at sunset. They are holding hands and walking towards the right. The woman is wearing a dark blue t-shirt, light-colored pants, a wide-brimmed hat, and sunglasses. The man is wearing a dark t-shirt, light-colored shorts, a baseball cap, and sunglasses. The ocean waves are breaking on the shore, and the sky is filled with soft, golden light from the setting sun. The couple's reflections are visible in the wet sand.

Interlocking Nail For
Cyrus Femur

Interlocking Nail For Cyrus Femur



Cyrus Femur Nail (L/R) Cat No. 187 / T187

Stainless Steel

Cat No.	Dia.	Length
187.L/R0934	09 MM	34 CM
187.L/R0936	09 MM	36 CM
187.L/R0938	09 MM	38 CM
187.L/R0940	09 MM	40 CM
187.L/R0942	09 MM	42 CM
187.L/R0944	09 MM	44 CM
187.L/R1034	10 MM	34 CM
187.L/R1036	10 MM	36 CM
187.L/R1038	10 MM	38 CM
187.L/R1040	10 MM	40 CM
187.L/R1042	10 MM	42 CM
187.L/R1044	10 MM	44 CM
187.L/R1134	11 MM	34 CM
187.L/R1136	11 MM	36 CM
187.L/R1138	11 MM	38 CM
187.L/R1140	11 MM	40 CM
187.L/R1142	11 MM	42 CM
187.L/R1144	11 MM	44 CM
187.L/R1234	12 MM	34 CM
187.L/R1236	12 MM	36 CM
187.L/R1238	12 MM	38 CM
187.L/R1240	12 MM	40 CM
187.L/R1242	12 MM	42 CM
187.L/R1244	12 MM	44 CM

Titanium

Cat No.	Dia.	Length
T187.L/R0934	09 MM	34 CM
T187.L/R0936	09 MM	36 CM
T187.L/R0938	09 MM	38 CM
T187.L/R0940	09 MM	40 CM
T187.L/R0942	09 MM	42 CM
T187.L/R0944	09 MM	44 CM
T187.L/R1034	10 MM	34 CM
T187.L/R1036	10 MM	36 CM
T187.L/R1038	10 MM	38 CM
T187.L/R1040	10 MM	40 CM
T187.L/R1042	10 MM	42 CM
T187.L/R1044	10 MM	44 CM
T187.L/R1134	11 MM	34 CM
T187.L/R1136	11 MM	36 CM
T187.L/R1138	11 MM	38 CM
T187.L/R1140	11 MM	40 CM
T187.L/R1142	11 MM	42 CM
T187.L/R1144	11 MM	44 CM
T187.L/R1234	12 MM	34 CM
T187.L/R1236	12 MM	36 CM
T187.L/R1238	12 MM	38 CM
T187.L/R1240	12 MM	40 CM
T187.L/R1242	12 MM	42 CM
T187.L/R1244	12 MM	44 CM



Cat. No. 179 / T179
Proximal Bolt 6.4 MM



Product Profile

Thread Dia.	:	6.4 MM
Core Dia.	:	4.8 MM
Head Dia.	:	8.0 MM
Screw Driver Use	:	5.0 MM

Cat. No. 151 / T151
I.L. Bolt 4.9 MM



Product Profile

Thread Dia.	:	4.9 MM
Core Dia.	:	4.3 MM
Head Dia.	:	8.0 MM
Screw Driver Use	:	4.5 MM

Cat. No. 182.C / T182.C
Cap For Cyrus Femur Nail

Size : 00 MM, 05 MM



Instrument Set For Cyrus Femur Interlocking Nail



Aiming Arm For Cyrus Femur Nail



Jig For Cyrus Femur Nail

Conical Bolt



Protection Sleeve (2 Nos.)



Drill Sleeve for 4.0 MM Drill Bit



Guide Wire Sleeve for 2.0 MM



Pointer Long



T - Handle (Wrench)



Spanner (Ordinary)
12 X 13, 18 X 13



Extractor Rod



Round Hammer



Impactor / Extractor Handle



Depth Gauge Extra Long



Bone Awl



Tissue Protector



Proximal Reamer 13 MM



Screw Driver 4.5 X 250 MM



C. C. Screw Driver 5.0 X 250 MM



Q.C. Drill Bits 6.4 MM



Drill Bits 4.0 MM X 10"



T Handle For
Q.C. Drill Bits



Interlocking Nail For
Intramedullary Supra Condylar
(I. M. S. C.)



Interlocking Nail For Intramedullary Supra Condylar (I. M. S. C.)



Interlocking I.M.S.C. Nail Cat No. 181 / T181

Stainless Steel

Cat No.	Dia.	Length
181.1015	10 MM	15 CM
181.1020	10 MM	20 CM
181.1025	10 MM	25 CM
181.1030	10 MM	30 CM
181.1115	11 MM	15 CM
181.1120	11 MM	20 CM
181.1125	11 MM	25 CM
181.1130	11 MM	30 CM
181.1215	12 MM	15 CM
181.1220	12 MM	20 CM
181.1225	12 MM	25 CM
181.1230	12 MM	30 CM
181.1315	13 MM	15 CM
181.1320	13 MM	20 CM
181.1325	13 MM	25 CM
181.1330	13 MM	30 CM

Titanium

Cat No.	Dia.	Length
T181.1015	10 MM	15 CM
T181.1020	10 MM	20 CM
T181.1025	10 MM	25 CM
T181.1030	10 MM	30 CM
T181.1115	11 MM	15 CM
T181.1120	11 MM	20 CM
T181.1125	11 MM	25 CM
T181.1130	11 MM	30 CM
T181.1215	12 MM	15 CM
T181.1220	12 MM	20 CM
T181.1225	12 MM	25 CM
T181.1230	12 MM	30 CM
T181.1315	13 MM	15 CM
T181.1320	13 MM	20 CM
T181.1325	13 MM	25 CM
T181.1330	13 MM	30 CM



Cat. No. 146 / T146

Cancellous 6.5 MM Screw (Fully Th.)



Product Profile

Thread Dia.	:	6.5 MM
Core Dia.	:	4.0 MM
Head Dia.	:	8.0 MM
Drill Bit Use	:	4.5 MM
Screw Driver Use	:	4.5 MM

Cat. No. 151 / T151

I.L. Bolt 4.9 MM



Product Profile

Thread Dia.	:	4.9 MM
Core Dia.	:	4.3 MM
Head Dia.	:	8.0 MM
Drill Bit Use	:	4.0 MM
Screw Driver Use	:	4.5 MM

Cat. No. 182 / T182

Cap For I.M.S.C. Nail

Size : 00 MM, 05 MM, 10 MM



Instrument Set For IMSC Interlocking Nail



Aiming Arm For I.M.S.C. Nail



Jig For I.M.S.C. Nail



Conical Bolt



Protection Sleeve



Pointer



Drill Sleeve for 4.0 MM Drill Bit



Drill Sleeve for 4.5 MM Drill Bit



T - Handle (Wrench)

Spanner (Ordinary)
12 X 13, 18 X 13



Extractor Rod



Round Hammer



Impactor / Extractor Handle



Depth Gauge Long



Bone Awl



Tissue Protector



Proximal Reamer 13 MM



Screw Driver 4.5 X 200 MM



Drill Bits 4.5 MM X 10"



Drill Bits 4.0 MM X 10"



A hiker wearing a red jacket and a backpack is seen from behind, walking up a narrow, rocky trail. The trail is heavily overgrown with lush green ferns and moss. The surrounding forest is dense with tall trees and more ferns, creating a vibrant green environment. The sky is visible at the top of the frame, showing a clear blue color.

Antegrade Femoral Nail
(A.F.N.)

Antegrade Femoral Nail



Antegrade Femoral Nail
Cat No. M196

Stainless Steel

Cat No.	Dia.	Length
M196.L/R09.34	09 MM	34 CM
M196.L/R09.36	09 MM	36 CM
M196.L/R09.38	09 MM	38 CM
M196.L/R09.40	09 MM	40 CM
M196.L/R09.42	09 MM	42 CM
M196.L/R09.44	09 MM	44 CM
M196.L/R10.34	10 MM	34 CM
M196.L/R10.36	10 MM	36 CM
M196.L/R10.38	10 MM	38 CM
M196.L/R10.40	10 MM	40 CM
M196.L/R10.42	10 MM	42 CM
M196.L/R10.44	10 MM	44 CM
M196.L/R11.34	11 MM	34 CM
M196.L/R11.36	11 MM	36 CM
M196.L/R11.38	11 MM	38 CM
M196.L/R11.40	11 MM	40 CM
M196.L/R11.42	11 MM	42 CM
M196.L/R11.44	11 MM	44 CM

Cat. No. M179

Proximal 6.4 MM Bolt



Product Profile

Thread Dia. : 6.4 MM
Core Dia. : 4.8 MM
Head Dia. : 8.0 MM
Screw Driver Use : 5.0 MM
Size : 24MM to 46MM

(in 2MM Variants)
50MM to

Cat. No. M151

I.L. Bolt 4.9 MM (Self Tepping)



Product Profile

Thread Dia. : 4.9 MM
Core Dia. : 4.3 MM
Head Dia. : 8.0 MM
Drill Bit Use : 4.0 MM
Screw Driver Use : 4.5 MM
Size : 60MM to 115MM

Variants)

(in 5MM

Instrument Set For Antegrade Femoral Nail



Main Jig For A.F.N



Aiming Arm



Conical Bolt



Screw Guide Device



AFN Bolt for Screw Guide Device (2 Nos.)



AFN Compression Screw



Slide Block For Hammer Unit



AFN Hammer Unit



AFN Removal of Nail



AFN Protection Sleeve 11 MM



Drill Sleeve 4.0 MM



Guide Wire Sleeve 2.5 MM



Long Pointer



Protection & Guide Sleeve 17 MM



T-Handle for Conical Bolt



Allen Key Driver for Conical Bolt

T-Handle With Q.C. for Reamer



Cannulated reamer for 6.4 MM Bolt (2.5 MM)



Drill Bits S.S. 12" X 4.0 MM



Screw Driver 5.0 MM Cannulated X 250 MM



Screw Driver for 4.5 MM X 250



Guide Pin Simple 16" X 2.5 MM



Guide Pin Threaded 16" x 2.5 MM



Guide Pin Simple 3.0 MM simple



Radiographic Scale



Femur Entry Awl CU. (T Type)



Spanner Ordinary S.S. 18 X 11



Proximal Reamer 13 MM



Depth Gauge Extra Long



Tissue Protector Simple



Measuring Device 2.5 MM




Flexible Reamer for Proximal 13.5 MM



Slotted Hammer







MEDI ORTHO LINK

Bringing Life Back To Your Bones

